

Chain of Custody

Standard Drinking Water Package

Customer Name: _____ Phone: _____

Address: _____ Email: _____

Sample Description: _____

Collected Date: _____ Collected Time: _____

NEL ID# _____

Preservation:

Bottle:

Analysis:

4°C

500 mL plastic

Alkalinity, Conductivity, pH,
F, Cl, NO₂, NO₃, PO₄, SO₄

HNO₃

250 mL plastic

As, Pb, Fe, Mn, Na, K, Ca, Mg, Hardness

Collected by: _____

Relinquished by: _____

Received by: _____

Date & Time: _____